



KEYSTONE CLEANROOM
PRODUCTS

STERILIZATION WRAP

*Providing an excellent bacterial and
microbial barrier for superior
performance post steam sterilization*


Ensure Maximum Compliance and Simplify your Operation

Keystone Cleanroom Products' line of sterilization wrappers are designed for steam sterilization. Depending on the material type selected, wrappers can be hydrophilic or hydrophobic and resist tearing and punctures.





PureClean™ 1421B
Tyvek® Sterilization
Wrap



BHD Sterilization Wrap



Heavy Duty SMS
Sequential
Sterilization Wraps



2-Ply SMS, One Step
Heavy Duty
Sterilization Wrap




2-Ply SMS, One-Step
Super Duty
Sterilization Wrap



Super Duty SMS
Sequential
Sterilization Wrap



1422A Tyvek®
Sterilization Wrap



GHD Wrapper



Custom solutions for your needs

Keystone Cleanroom Products offers tailored sterilization wrapping solutions designed to meet the specific needs of your operation. Choose from a variety of materials, including Tyvek® (1422A or 1421B), SMS, GHD, or BHD.

Temperature Requirements for Steam Sterilization By Material

When choosing the right material for your operation, it is important to understand the temperature limits that each material can withstand. See the chart below to learn more.

Material Type	Temperature Limitations
1421B/1422A Tyvek®	Steam sterilization cycles not to exceed 250°F to 260°F [121°C to 127°C] at 30 psi, for 30 minutes.
SMS	Steam sterilization cycles not to exceed 275°F (135°C)
BHD/GHD	Steam sterilization cycles not to exceed 275°F (135°C)



SINCE 1925

KEYSTONE

SAFETY ✧ CLEANROOM PRODUCTS

Contact Us:

info@keystonecrp.com

+1 856-663-5740