


KEYSTONE CLEANROOM  
PRODUCTS



# 1421B TYVEK® PRODUCT CATALOG

*Designed with Annex 1 Compliance  
in mind*

# A line of products designed with Annex 1 in mind

Partner with us to achieve complete Annex 1 compliance through materials and consultation that fully adhere to all regulatory standards. We also deliver tailored solutions, offering custom sterilization wrapping designed to meet your operation's requirements.

Additionally, we provide ongoing support, including validation assistance and post-project consultation, to ensure continued compliance and success.

---





## Autoclave Bags

- Heat Seal, Self- Seal, and Drawstring Closures Available
- Contact us for sizing details



## Autoclavable Stopper Bowl Covers

- Contact us for sizing details



## Autoclavable Muff Covers

- Contact us for sizing details



## Autoclavable Continuous Form Tubing

- Contact us for sizing details



## Autoclavable Wrappers

- Contact us for sizing details



## Autoclavable Filling Needle Bags

- Contact us for sizing details



SINCE 1925

**KEYSTONE**

SAFETY ♦ CLEANROOM PRODUCTS




## PureClean™ 1421B DuPont™ Tyvek®

Keystone is proud to partner with DuPont™ to offer the new PureClean™ 1421B Tyvek® material as a potential replacement for the legacy 1422A type. PureClean™ Tyvek® is essentially the same as the original but without the antistat agent. This advanced material is designed for pharmaceutical manufacturing customers seeking the benefits of the next generation of Tyvek®.

SINCE 1925

**KEYSTONE**

SAFETY ↔ CLEANROOM PRODUCTS



## Clear Film Side for Visual Inspection

Keystone Cleanroom Products'

PureClean™ Self-Sealing 1421B Tyvek®  
Autoclave bag with Clear Film Side and  
Steam Indicator has a clear film side that  
allows for easy visual inspection without  
breaking the sterile barrier, as required by  
Annex 1.



## Hermetically Sealed

Keystone Cleanroom Products provides autoclave bags with both heat-seal and self-seal closures, eliminating the challenges posed by the drawstring closures found in larger Tyvek bags. These drawstring closures often make it difficult to fully dry the folds and creases of the twisted "gooseneck" design.



## "Right-sizing" Bags Used in Production

If your operation isn't ready to switch to the new 1421B PureClean Tyvek, Keystone Cleanroom Products can offer some solutions while still using the legacy 1422A Tyvek. Using the correct size self-sealing 1422A Tyvek bag for steam sterilization, instead of an oversized drawstring bag, eliminates the excess hydrophobic material; and reduces the chance of any trace of the water soluble anti-stat agent residue being left behind.




SINCE 1925

**KEYSTONE**

SAFETY ⇄ CLEANROOM PRODUCTS

# Optimizing Compliance: Keystone's Innovative Solutions

We have the ability to design customer specific sterilization wrapping systems and manufacture our products in Pennsauken, NJ USA. We operate in certified ISO6 and ISO8 cleanrooms, with only the most advanced manufacturing equipment, materials, and processes.





SINCE 1925

**KEYSTONE**

---

SAFETY ✦ CLEANROOM PRODUCTS

**Your Journey to Annex 1  
Compliance in the Steam  
Sterilization of  
Pharmaceutical Change  
Parts Begins Here.**

Contact Us:  
[info@keystonecrp.com](mailto:info@keystonecrp.com)  
+1 856-663-5740